

# NONWOVEN DISPOSABLE GOWNS

**AAMI -LEVEL 2**

**OPTIONS**

STERIL / NON STERIL

**REF: SSU-902**

Item name	Nonwoven Disposable Gowns
Style	SS Surgical Gown with knit cuff
Material	SS 40 gsm (Spunbond+Spunbond)
Size	M-L-XL-XXL
Colour	Blue
Certificate	CE,FDA

## Level 2

Used in LOW risk situations

Provides a barrier to larger amounts of fluid penetration through splatter and some fluid exposure through soaking

Two tests are conducted to assess barrier protection performance:


Water impacting the surface of the gown material

Blood draw from a vein, Suturing, Intensive care unit, Pathology lab



# NONWOVEN DISPOSABLE ISOLATION GOWNS

**AAMI -LEVEL 4**



**REF: SSU-904**

Item name	Nonwoven Disposable Isolation Gown
Style	SMS Isolation Gown with knit cuff
Material	Material: SMS 40 gsm + PE Breathable 19 gsm Laminated Film
Size	M-L-XL-XXL
Colour	Blue
Certificate	CE,FDA




**Level 4**

Used in HIGH risk situations

Prevents all fluid penetration for up to 1 hour

May prevent VIRUS penetration for up to 1 hour  
In addition to the other tests conducted under levels 1-3, barrier level performance is tested with a simulated blood containing a virus. If no virus is found at the end of the test, the gown passes.

Pathogen resistance,  
Infectious diseases (non-airborne),  
Large amounts of fluid exposure over long periods



## PRESSURE TESTS

